





Reduce Risk with Right Sign On-Site Coronavirus Testing








Requires Only a Drop of Blood

The lgG/lgM Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of lgG and lgM antibodies to SARSCoV-2 in human whole blood, serum, or plasma. The Rapid Test is an effective means to detect and inform if employees or workers have had COVID-19 or if they are free from the virus.


Process



PRODUCT # FCPPE-RAPID

▲ COVID-19 IgG/IgM Rapid Test is for professional use. It is NOT intended for self-testing or home use.

Accurate Results



RAPID & PCR TEST KITS IN STOCK & READY TO SHIP TODAY!

MINIMIZE RISK WITH CUSTOM ON-SITE TEST PROGRAMS









How COVID-19 Tests Work Together

Dedicated testing & monitoring is the most effective way to minimize risk of a COVID-19 outbreak for yourself & your employees

	RAPID TEST	PCR TEST
How it Works	Detects the body's immune response to the virus, called antibodies (IgM, IgG)	Detects the presence of the virus' genetic material (RNA)
Type of Sample Needed	Blood sample (Finger Prick or Blood Draw)	Nasopharngeal (nose/throat) swab
Time for Results	10 - 15 minutes	24-72 hours after receipt from lab
When to Use	Screening populations to determine infection risk, and asymptomatic people	If exhibiting symptoms or positive on an Antibody test
What a Positive Result Means	Patient has/had the virus in the past 3 months	Patient is infected with live virus

First Call PPE works with top clinicians, employers, manufacturers, labs, & pharmacists to develop and execute practical, scalable, & sustainable testing solutions.